

Bio Light 3

The Light of Life



THE REALITY

Research has shown that all living cells of plants, animals and human beings emit biophotons – electromagnetic waves of light within biological systems – which cannot be seen by the naked eye but can be measured using special equipment. Healthy cells, sick cells, tumour cells, and infectious viruses and bacteria all have different light signals, which come together to form an internal system of cellular communication. A healthy cell radiates coherent light; whereas an unhealthy cell radiates chaotic light.

OUR DISCOVERY

The Amezcua Bio Light 3 is a breakthrough in biophoton light therapy – a holistic approach to wellness that covers the physical, emotional, mental and spiritual elements of well-being. Now, no matter where you are, you can get the benefits of light in a safe environment.

THE RESULTS

Since the body and mind are connected through our bio-energy, not only can this therapy improve physical conditions, it also helps in achieving a more harmonised state of mental and emotional health.

THE BENEFITS OF EACH COLOUR

Red

- Elevates mood
- Activates hormones
- Promotes cell regeneration

Blue

- Anti-inflammatory qualities
- Increases metabolism and builds vitality
- Calms nerves
- Regulates sleep pattern

White

- Purification of space
- Balance and harmony
- Raises the vibrations of one's consciousness & the body

Green

- Aids in detoxification
- Relieves stress
- Boosts the immune system

Near Infrared (Always on)

- Aids in wound healing
- Improves blood circulation
- Eases pain
- Helps with skin repair



HOW TO USE

1. Open the front lid and place the Amezcua Bio Disc 3 inside the Amezcua Bio Light 3 compartment.
2. Touch and hold the ON sensor for 5-7 seconds and shine the Bio Light 3 on the desired part of your body. Avoid using the Amezcua Bio Light 3 on the eyes, metal implants fixed to joints or bones, over the brain (in case of malignant tumors), head area for individuals with epilepsy conditions, on the foetus area (if pregnant) and over the heart (in case of any heart diseases).
3. Choose the light color mode according to your need. (refer to "THE RESULTS")
4. Use this for 3-4 sessions a day, 15 minutes per session.

GET YOUR BIO LIGHT 3 TODAY AND AWAKEN THE LIGHT WITHIN!

Disclaimer:

The information in this above does not replace medical advice. Advice on treatment or care of an individual patient should be obtained through consultation with a physician or trained health care practitioner who has examined that patient or is familiar with that patient's medical history.

Please visit www.qnet.net for more info and references.